



FARM SAMPLING

PROCESS SAMPLING

PERFORMANCE ANALYSIS

## 4 oz High Profile Vial

One-Piece Flip-Top Vial - Specific for aseptic sampling



### Applications

#### Ideally shaped for

- Field collection of farm milk and related dairy products
- Veterinary specimen collection, infection control of mastitis
- Environmental sampling containers

#### Options

- Can be supplied with Bronopol and / or BARCODE label
- Standard Colour is natural white – other colours on request with minimum production run requirements

### Benefits

#### Aseptic single-use vial with flip top construction

- Manufacture in clean room
- Patented « in mould closing » process
- Perfectly tight and aseptic
- Translucent in colour for visual inspection
- Pure Polypropylene suitable for food contact
- 100% Recyclable after use

#### User-friendly hinged vial concept

- One hand opening and closing for perfect cleanliness and ease of use
- No risk of contamination due to polluted stoppers
- Allows the use of gloves in cold environment

#### Safety and Quality

- Alimentary : Complies with CE Norm N° 1935/2004
- Traceability of batches of production : a label on each box includes batch number and production date

### Characteristics

Internal Volume	130,0 ml
Vial Overall Height	90,8 mm
Bottom Body Diameter	45,2 mm
Top Body Diameter	46,7 mm
Top Lid Diameter	49,8 mm
Body Height	89,3 mm
Inside Top Diameter	44,6 mm
Inside Height	88,3 mm
Level Mark Line	95,0 ml
Vial Weight	17,6 grams

PRODUCTION MOULD # 800

ART. # 03 04 7090 00HP CN

In colour Clear - 250 Pieces per box - 30 Boxes - 7.500 Pieces per Pallet