

TECHNICAL DATA

DESCRIPTION

- Filter card for ADMI test composed of a non-woven filter
- Material: 90% polyester 10% cellulose
- Binder: yes

APPLICATIONS

- Determination of burnt particles in milk powders (according to ADMI standard)
- Determination of sediments in food products

TECHNICAL SPECIFICATIONS

CHARACTERISTICS	VALUES	UNITS	STANDARDS
Weight	110	g/m ²	DIN 53104
Thickness	0,7	mm	
Porosity	50	µm	
Air permeability	2200	l/m ² /s	
Water flow	510	l/m ² /s	
M.D elongation resistance	113	N/5 cm	
C.D elongation resistance	77	N/5 cm	

COMMERCIAL DATA

SIZES	REFERENCES	PACKING
57x155 mm	09.0110	500
45x80 mm	06.8050	500