

|                     |                                          |        |                  |                  |                                                                                     |
|---------------------|------------------------------------------|--------|------------------|------------------|-------------------------------------------------------------------------------------|
| PDS No. 7222xx      | PRODUCT DATA SHEET                       |        |                  |                  | Page 1 of 1                                                                         |
| Revision 07         | Reaction Tubes, 2 ml, Free-Standing Base |        |                  |                  |  |
|                     | Greiner Item-No. 7222xx                  |        |                  |                  |                                                                                     |
| Valid for Item-No.: | 722201                                   | 722251 | 722261 (sterile) | 722263 (sterile) |                                                                                     |

| 1.  | Description / Specification |                                                                                                                                                                                                                                                                                                                                                                  |
|-----|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1.1 | Description                 | Reaction tube, conical, 2 ml, with graduation, writing area and free-standing base<br>722201: without screw cap, non-sterile<br>722251: without screw cap, sterile<br>722261, -263: with screw cap, sterile<br><u>Appropriate screw caps:</u><br>Cat.-No. 366380, translucent<br>366381, white<br>366383, red<br>366384, blue<br>366385, green<br>366386, yellow |
| 1.2 | Dimensions                  | See Customer Drawing                                                                                                                                                                                                                                                                                                                                             |
| 1.3 | Volume                      | 2.0 ml                                                                                                                                                                                                                                                                                                                                                           |
| 1.4 | Material / Resin            | <u>Tube:</u> PP (Medical Polypropylene), free of heavy metal<br><u>Cap:</u> HDPE (High Density Polyethylene, medical grade), free of heavy metal<br><u>"O"-Ring (Cap):</u> thermoplastic synthetic rubber, free of heavy metal                                                                                                                                   |
| 1.5 | Colour                      | <u>Tube:</u> translucent<br>722261: <u>cap:</u> translucent<br>722263: <u>cap:</u> red                                                                                                                                                                                                                                                                           |
| 1.6 | Sterilization               | 722201: no<br>722251, -261, -263: SAL10 <sup>-3</sup>                                                                                                                                                                                                                                                                                                            |
| 1.7 | Quality Control             | <u>Product-Control:</u> testing of attributive and variable characteristics in accordance with the valid specification                                                                                                                                                                                                                                           |
| 1.8 | Intended Use                | General laboratory product for the processing, centrifugation and storage of samples to be used by qualified personnel in a laboratory environment.                                                                                                                                                                                                              |
| 1.9 | Other Information           | For single use only                                                                                                                                                                                                                                                                                                                                              |

| 2.  | Features                 |                                                                                                                                                                    |
|-----|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 2.1 | Basic features           | Free of detectable DNase/RNase and human DNA, non-pyrogenic                                                                                                        |
| 2.2 | Temperature range        | For application: -80°C to +121°C                                                                                                                                   |
| 2.3 | Autoclavability          | Yes                                                                                                                                                                |
| 2.4 | Centrifugation, max. RCF | 22.000 x g: fixed-angle rotor                                                                                                                                      |
| 2.5 | Chemical Resistance      | See homepage:<br><a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a> |
| 2.6 | Shelf life               | 722201: n/a<br>722251, -261, -263: 5 years                                                                                                                         |
| 2.7 | Other Information        | Store in dry place, temperature range: -10°C to 50°C                                                                                                               |

| 3.  | Packaging         | 722201, -251                                                       | 722261 | 722263 |
|-----|-------------------|--------------------------------------------------------------------|--------|--------|
| 3.1 | Pieces / Bag      | 500                                                                | 500    | 50     |
| 3.2 | Pieces / Box      | 5.000                                                              | 3.000  | 3.000  |
| 3.3 | Lot-No.           | YY MM XXXXXX (year, manufacturing facility, internal order number) |        |        |
| 3.4 | Other Information | -                                                                  |        |        |

| 4. | Other Information |
|----|-------------------|
|    | -                 |

Data Sheet subject to change without notice!

| Prior Issue        | Drawn                    | Approved                 | Released                 | <b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved. |
|--------------------|--------------------------|--------------------------|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Revision 06        | Date<br>1 September 2021 | Date<br>24 November 2021 | Date<br>25 November 2021 |                                                                                                                                                                                                                                                                                                       |
| Date<br>04.03.2019 | Name<br>S. Kaelberer     | Name<br>Dr. T. Schreiber | Name<br>A. Illig         |                                                                                                                                                                                                                                                                                                       |

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.