



Item Number	Title	Quantity
1	WST1 Shelf	1
2	650x50x25 Box Section	2
3	610x50x25 Box Section	2
4	RVN-0005 M8 Riv Nut	12
5	3 inch castor - 367H00G080XZPF	4

Revision	Chg Doc	Date	Tolerances Unless Otherwise Specified Fabrication ± 0.30 Machining ± 0.05 Angles $\pm 0.50^\circ$ Dimensions in mm Remove all burrs and sharp edges Third Angle Projection	Material Mild Steel Finish Powder Coated Size A3 Computer Generated Drawing Do Not Revise Manually Do Not Scale	Drawn CS Date 11/03/09 Shr 1 of 1 Approved Date This document and the information contained in it must not be disclosed and must not be reproduced in whole or in part without the written authorisation of Leec Limited	LEEC LEEC LIMITED PRIVATE ROAD No 7 COLWICK INDUSTRIAL ESTATE NOTTINGHAM, NG4 2AJ		
						Description WST1 Wheeled Platform		
						Client		
						Drg No.	Rev A	Scale