



Be Right™



## Ammonium cuvette test 0.5-5.0 mg/L NH<sub>4</sub>-N, 25 tests

Product #: LCK505

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

---

### Specifications

According to standard:	ISO 7150-1, DIN 38406 E5-1, UNI 11669:2017
Instrument:	DR1900, DR3900, DR6000
Measuring range:	0.5 - 5.0 mg/L NH <sub>4</sub> -N
Measuring range (2):	0.65 - 6.4 mg/L NH <sub>4</sub>
Method:	Indophenol Blue
Number of tests:	25
Parameter:	Ammonium
Platform:	LCK
Shelf life:	17 months from production date
Storage conditions:	2 - 8 °C (keep refrigerated)