



# CSP Dairy™ 2 oz Solid-Top Vial

Single Use Clean Flip-Top Vial for Dairy,  
Food and Veterinary Sampling

## Applications

### Fields of use

- Field collection of farm milk and related dairy products
- Quality control sample storage containers and food packaging vials
- Veterinary specimen collection, infection control of mastitis
- Environmental sampling containers

### Options

- Can be supplied with bronopol
  - Ready to use upon request
  - Fits perfectly into common transport racks and trays for DHIA
- Can be supplied with label (barcode...)
- Standard vial is colourless - colours upon request

## Benefits

### Flip-Top Design

- Consumer-friendly concept (easy open/close)
- Patented “in mould closing” process
- Perfectly tight and clean by process
- Translucent in colour for visual inspection
- Pure Polypropylene suitable for food contact
- 100% Recyclable after use

### User-Friendly Hinged Vial Concept

- One-hand opening and closing for perfect cleanliness and ease of use
- Allows the use of gloves in cold environment

### Regulatory and Quality

- ISO-9001 certified company
- Complies with Regulation (EC) No 1935/2004

## Characteristics

Please note that these values are average measurements  
Tolerance +/- 0.5 mm (will vary based on product)

Internal volume	49 ml
Filling level line	33 ml
Vial overall height	84,7 mm
Body height	82,5 mm
Inside height	81,4 mm
Top lid diameter	35,1 mm
Top body diameter	30,9 mm
Inside top diameter	29,1 mm
Bottom body diameter	28,9 mm
Vial weight	10,5 grams

Production mould # 728

## Packaging

ART. # 03 03 1084 00LH CN

- Clear in color
- 600 Pieces per box
- 30 Boxes
- 18.000 Pieces per pallet

