

BagPipet®

Fixed-volume pipet & straw tips

Ref 251 091 (BagPipet®)



Standards:



Manufactured under:



BagPipet® is an autoclavable pipet with fixed volume (0.1 mL / 0.9 mL / 1 mL). Used with BagTips® a single use rigid straw tip, available in 2 sizes (Regular 19 cm / Jumbo 24 cm), BagPipet® intakes any product, even very viscous products: milk, cream, yoghurt, cosmetics...The long BagTips® reaches the bottom of any bag.

TECHNICAL SPECIFICATIONS

- Fixed volume pipet
- Ejection of disposable straw
- Light and ergonomic
- Autoclavable
- Easy calibration
- 3 settings: 0.1 mL / 0.9 mL / 1 mL
- Calibrated according to **the GLP "Good Laboratory Practice"**
- Also available: volume push button
 - 100 µL (Ref 251 001)
 - 900 µL (Ref 251 009)
 - 1000 µL (Ref 251 010)
 - 100 à 1000 µL (Ref 251 080) Variable
- Dimensions : 6 x 18,5 x 3 cm, net weight: 0,149 kg
- Box (w x d x h): 30 x 10 x 6 cm, gross weight: 250 g
- 1 year warranty
- Availability of spare parts: 10 years
- In compliance with: ISO 7218, ISO 6887, ISO 4833-2 and FDA BAM (Bacteriological Analytical Manual)
- Designed and made in France

DELIVERED WITH

- 1 user's manual
- 1 pouch BagTips® Regular
- 1 holder for key
- 1 volume button 1000 µL (mounted on the pipet)
- 2 volume buttons (900 µL, 100 µL)
- 2 O-Rings
- 1 Allen wrench
- Oil

ACCESSORIES

- Ref. 252 019 BagTips® Regular : Straw tip (5 to 12 kGy gamma ray treated) for BagPipet.
Nominal volume: 1.1 mL - Length: 190 mm – Quantity per box: 1000
Box (w x d x h): 27,5 x 23 x 11 cm, gross weight: 1,2 kg
- Ref. 252 024 BagTips® Jumbo : Straw tip (5 to 12 kGy gamma ray treated) for BagPipet.
Nominal volume: 2,2 mL – Length: 240 mm - Quantity per box: 1000
Box (w x d x h): 27,5 x 23 x 11 cm, gross weight: 1,5 kg