

Dilucup Elegance BPW (Buffered Peptone Water)

Use and description: Ready-to-use plastic cups, each contains 9 mL medium to be used for diluting microbiological samples together with a special shaker, Dilushaker or Diligent Shaker. By adding a sample of 1 mL to the first cup a dilution of 1:10 will be achieved. 1 mL from this cup is then transferred to the next to get a dilution of 1:100. The procedure is continued until the required dilution is reached.

Composition per litre:

Peptones	10,00 g
NaCl	5,00 g
Buffers	3,56 g

pH: 7,0 +/-0,2 at 25°C

Volume per cup: 9 mL +/- 2 %

Sterilization: E-Beam

Material of the cups: PET

Expiry date: When unopened 30 months after production.

Storage: Unopened: +2 to +30°C. +2 to +8°C when the box is open.

Packaging: The cups are delivered in blisters of 7 pre-cut rows of 3 cups (40DEBPW3) or 6 cups (40DEBPW6). Each box contains 672 cups distributed in 2 bags of 336 cups. Weight and volume: 8,1 kg/box - 0,032 m³/box

DILUCUP ELEGANCE BPW 2 X 16 X 21 CUPS

ARTICLE NUMBER: 40DEBPW3

DILUCUP ELEGANCE BPW 2 X 8 X 42 CUPS

ARTICLE NUMBER: 40DEBPW6

