

## CERTIFICATE OF STERILITY AND QUALITY ASSURANCE

PRODUCT	LOT NUMBER	EXPIRY DATE	CERTIFICATION DATE
M452022			

### Environmental-Large -Sponge swab-Neutralising Buffer.

PARAMETER	SPECIFICATION	RESULT
Sterilisation Method	Irradiation	PASSED
Sterility Test	Absence of aerobes and anaerobes after immersion in tryptic soy broth for 10 days at 30°C	PASSED
Performance Testing	<b>Recovery within specification of following organisms:</b> <i>Escherichia coli</i> NCTC 9001 (24 hr) <i>Salmonella typhimurium</i> NCTC 0074 (24 hr) <i>Staphylococcus aureus</i> NCTC 8532 (24 hr) <i>Listeria monocytogenes</i> NCTC 5214 (24 hr)	PASSED

Product has been manufactured in accordance with ISO 9001:2015

Performance testing conducted in accordance with: Internal Procedures

For and on behalf of Medical Wire & Equipment



Quality Assurance