


Biodegradable Nitrile Powder Free Gloves (100 per Pack)



Size & Dimensions	Palm Width (mm)	Length (mm)
S	<85	240
M	<95	240
L	<105	240

- Finger Textured
- Not made with natural rubber latex
- Ambidextrous
- Patient examination glove
- Single use only



Colour	Violet Blue
Powder/Powder-Free?	Powder-Free
Thickness	3.5g
Beaded Cuff	Yes
Material	Nitrile
Carton Qty	10 Boxes/Carton

PHYSICAL PROPERTIES

	Before Aging	After Aging
Tensile Strength	≥18 MPa	≥16 MPa
Elongation	≥500%	≥400%
Strength at Break	≥6N	≥6N