



Reference A04050S

One pair of plastic over-shoes + one pair of gloves



Standard ISO 13307

Microbiology of food and animal feed
Primary production stage Sampling techniques



European Union
regulations

Use:

Used for Salmonella testing in poultry farming,
for sampling surfaces under the boot.

To be used with sterile SteriSox® references:
A04130, A04137, A04139, A04142 and A04149

Features:

- 1 pair of plastic over-shoes
 - Universal size
 - With elastic band at top
- 1 pair of blue gloves
 - Vinyl quality
 - Size 9 (L)



Reference	Custom code	Packaging	Cardboard box size
A04050S	3926200000	100 units per box packed unitary in a zip bag size 18 x 25 cm	60 x 40 x 25 cm

Sterility and traceability:

Sterilization	Sterilization certificate	Shelf life	Traceability *
Sterilization by gamma rays Calibrated dose minimum 20 kiloGrays	in each box	24 months	On each product/box



Before



After

* Printed information: reference, designation, batch number, diluent code, production and expiry date

04/01/2024