



## Reference A04730

**SteriSox® non-woven, 1 pair, dry in a zip sampling bag**



### ISO 13307 Standard

Microbiology of food for human and animal consumption. Primary production stage. Sampling techniques.



**European regulation**

### Use:

Sterile surface control device for collecting dust and feces from the floor.

Microbiological monitoring of Salmonella during rearing:

- E.g., birds subject to European regulations: *Gallus gallus*, turkeys
- E.g., birds not subject to European regulations (duck, guinea fowl...) and fish, pig...

### Features:

- **SteriSox® non-woven**
- 1 pair of boots swabs
- With elastic band
- **Diluent**
- No impregnation (dry)
- **Zip sampling bag**
- Tearable with zipper, 17.5 × 27.5 cm
- Identification area



2 SteriSox® non-woven



Product code	Packaging
A04730	100 units per box

### Sterility and traceability:

Sterilization	Sterilization certificate	Validity	Traceability *
by gamma rays Calibrated dose minimum 20 kiloGrays	On each carton	2 years from the sterilization date	Printed on each unit/box



Before



After

\* Printed information: reference, designation, batch number, diluent code, production and expiry date

03/12/2025

