



SteriSox® non-woven



Reference A04736

Sterisox® non-woven, 1 pair, tryptone salt in zip sampling bag, with 2 gloves and 2 plastic overshoes



ISO 13307 Standard

Microbiology of food for human and animal consumption. Primary production stage. Sampling techniques.



European regulation

Use:

Sterile surface control device for collecting dust and feces from the floor.

Microbiological monitoring of Salmonella during rearing:

- E.g., birds subject to European regulations: *Gallus gallus*, turkeys
- E.g., birds not subject to European regulations (duck, guinea fowl...) and fish, pig...

Features:

- **Sterisox® non-woven**
 - 1 pair of boots covers
 - Retention elastic
- **Diluent**
 - Tryptone salt
 - Calibrated dose: 15 ml
- **Zip sampling bag**
 - Tearable with zipper, 17,5 x 27,5 cm
 - Identification area
- **Sampling accessories**
 - 1 pair of plastic overshoes: universal size, with elastic band at top
 - 1 pair of gloves: blue, polyethylene, ambidextrous, latex-free, powder-free, universal size, packed in a small bag 13 x 10 cm
 - Packed together in an outer zip bag with the sterile sampling device SteriSox®



Product code	Packaging
A04736	65 units per box

Sterility and traceability:

Sterilization	Sterilization certificate	Validity	Traceability *
by gamma rays Calibrated dose minimum 20 kiloGrays	On each box	1 years from the sterilization date	Printed on each unit/box



* Printed information: reference, designation, batch number, diluent code, production and expiry date

03/12/2025

