



Reference A04739

SteriSox® non-woven, 1 pair, buffered peptone water at 10% neutralizing agent in zip sampling bag



ISO 13307 Standard

Food Food microbiology for human and animal consumption. Primary production stage. Sampling techniques



European regulation

Use:

Sterile device for microbiological control of Salmonella on surfaces. Allows verification of cleaning and disinfection effectiveness, even in the presence of disinfectant residues.

Features:

- **SteriSox® non-woven**
 - 1 pair of boot swabs
 - With elastic band
- **Diluent**
 - Buffered peptone water at 10% neutralizing agent
 - Calibrated dose: 15 ml
- **Zip sampling bag**
 - Tearable with zip, 17,5 × 27,5 cm
 - Identification area



Product code	Packaging
A04739	90 units per box

Sterility and traceability:

Sterilization	Sterilization certificate	Validity	Traceability *
by gamma rays Calibrated dose minimum 20 kiloGrays	On each box	1 year from the sterilization date	Printed on each unit/box



Before



After

* Printed information: reference, designation, batch number, diluent code, production and expiry date

03/12/2025

