



Reference A04742

SteriSox® non-woven, 2 pair, tryptone salt in zip sampling bag



ISO 13307 Standard

Food microbiology for human and animal consumption.
Primary production stage.
Sampling techniques.



European Regulation

Use:

Sterile surface control device for collecting dust and feces from the floor.

Microbiological monitoring of Salmonella during rearing:

- E.g., birds subject to European regulations: *Gallus gallus*, turkeys
- E.g., birds not subject to European regulations (duck, guinea fowl...) and fish, pigs

Features:

- **SteriSox® non-woven**
 - 2 pairs of boot swabs
 - With elastic band
- **Diluent**
 - Tryptone salt
 - Calibrated dose: 30 ml
- **Sampling Bag**
 - Tearable with zip, 17,5 × 27,5 cm
 - Identification area



Product code	Packaging
A04742	60 units per box

Sterility and traceability:

Sterilization	Sterilization certificate	Validity	Traceability *
by gamma rays Calibrated dose minimum 20 kiloGrays	On each carton	1 year from the sterilization date	Printed on each unit/box



Before



After

* Printed information: reference, designation, batch number, diluent code, production and expiry date

03/12/2025

