



SteriSox® non-woven

## Reference A04745

**SterisoX® non-woven, 5 pairs, tryptone salt in 2 zip sampling bags, with 2 gloves and 2 plastic overshoes**



**European regulation**



**ISO 13307 Standard**

Food microbiology for human and animal consumption. Primary production stage. Sampling techniques

### Use:

Sterile surface control device for collecting dust and feces from the floor.

Microbiological monitoring of Salmonella during rearing:

- E.g., birds subject to European regulations: *Gallus gallus*, turkeys
- E.g., birds not subject to European regulations (duck, guinea fowl...) and fish, pigs

### Features:

- **SterisoX® non-woven**
  - 5 pairs of boot swabs
  - With elastic band
- **Diluent**
  - Tryptone salt
  - Calibrated dose: 70 ml (2 x 35 ml)
- **2 zip sampling bags** (5 SterisoX® non-woven per bag)
  - Tearable with zip, 17,5 x 27,5 cm
  - Identification area
- **Sampling accessories**
  - 1 pair of plastic over-shoes: universal size, with elastic band at top
  - 1 pair of gloves: blue, polyethylene, ambidextrous, latex-free, powder-free, universal size, packed in a small bag 13 x 10 cm
  - Packed together in an outer zip bag with the sterile sampling device SterisoX®



Product code	Packaging
A04745	20 units per box

### Sterility and traceability:

Sterilization	Sterilization certificate	Validity	Traceability *
by gamma rays Calibrated dose minimum 20 kiloGrays	On each box	1 year from the sterilization date	Printed on each unit/box



\* Printed information: reference, designation, batch number, diluent code, production and expiry date

03/12/2025

