



Reference A04748

SterisoX® non-woven, 2 pairs, tryptone salt in a zip sampling bag, with 2 gloves and 2 plastic overshoes



ISO 13307 Standard

Food microbiology for human and animal consumption.
Primary production stage.
Sampling techniques



European regulation

Use:


Sterile surface control device for collecting dust and feces from the floor.

Microbiological monitoring of Salmonella during rearing:

- E.g., birds subject to European regulations: *Gallus gallus*, turkeys
- E.g., birds not subject to European regulations (duck, guinea fowl...) and fish, pigs

Features:

- **SterisoX® non-woven**
 - 2 pairs of boot swabs
 - With elastic band
- **Diluent**
 - Tryptone salt
 - Calibrated dose: 30 ml
- **Zip sampling bag**
 - Tear-open with zipper, 17,5 x 27,5 cm
 - Identification area
- **Sampling accessories**
 - 1 pair of plastic over-shoes: universal size, with elastic band at top
 - 1 pair of gloves: blue, polyethylene, ambidextrous, latex-free, powder-free, universal size, packed in a small bag 13 x 10 cm
 - Packed together in an outer zip bag with the sterile sampling device SteriSox®



Product code	Packaging
A04748	50 units per box

Sterility and traceability:

Sterilization	Sterilization certificate	Validity	Batch Number
Gamma ionization Minimum 20 kilograys	On each box	1 year from the sterilization date	Printed on each unit/box



Before



After

* Printed information: reference, designation, batch number, diluent code, production and expiry date

03/12/2025

