

Powerhold Diamond Nitrile Grip Powder-Free Gloves Orange

Box Quantity: 100pcs
(1000pcs per outer case)



KEY FEATURES

- Diamond pattern for superior grip
- Over double thickness compared to standard nitrile
- PPE Category III
- Beaded cuffs
- Powder-free and latex-free
- Food-Safe EN 1186
- Popular with automotive and industrial
- Provides comfort and protection from hazardous chemicals
- Unique orange colour for visibility



TECHNICAL INFORMATION

Material	Nitrile	Colour	Orange
Available Sizes	S, M, L, XL	Sterile	Non-Sterile
Glove Length (mm)	≥240	Freedom from Holes	AQL 1.5
Powder Content	Powder-Free	Surface	Diamond grip texture
Finger Thickness (mm)	0.25 ± 0.03	Palm Thickness (mm)	≥0.22
Force at Break	≥6N before and after aging	Cuff	Beaded Cuff
Elongation	≥500% before aging		
Product Code	Size & Dimensions	Palm Width (mm)	Length (mm)
INC501854	Small	<84	240
INC501857	Medium	<94	240
INC501860	Large	<105	240
INC501863	Extra Large	<114	240

Packaging & Storage

Packaging	Case 10 boxes of 100
Box Size (mm)	230 x 125 x 90
Carton Dims (mm)	465 x 255 x 265
Shelf Life	Expiration date is 5 years after production date
Storage Instructions	Keep out of direct sunlight, store in a cool & dry place. Keep away from sources of ozone or ignition.
Accreditations	EN ISO 14971:2019; EN ISO 15223-1: 2016; EN 1041:2008+A1:2013; ISO 10993-1: 2018; EN ISO 10993-5: 2009 EN ISO 10993-10: 2013; EN455-1:2020; EN455-2:2015; EN455-3: 2015; EN455-4:2009; EN ISO 21420:2020 EU 2017/745 Class I; EU 2016/425