

## Vertex Nitrile Powder-Free Gloves Green Long Cuff

Box Quantity: 100pcs (1000pcs per outer case)





## **KEY FEATURES**

- Strong and durable glove to prevent tearing in heavy duty work
- Latex free and Powder Free for sensitive skin
- Comply with PPE CAT III regulation and EN455 rated
- Suitable for use in a wide range of industries including automotive, life sciences, food processing and chemical handling.
- · Micro-texture over full surface



























TECHNICAL INFORMATION			
Material	Nitrile	Colour	Green
Available Sizes	S, M, L, XL	Sterile	Non-Sterile
Glove Length (mm)	≥290	Freedom from Holes	AQL 1.5
Powder Content	Powder-Free	Surface	Micro-textured on whole surface
Finger Thickness (mm)	0.18 ± 0.03	Palm Thickness (mm)	≥0.14
Product Code	Size & Dimensions	Palm Width (mm)	Length (mm)
INC502055	Small	<85	290
INC502058	Medium	<96	290
INC502061	Large	<107	290
INC502064	Extra Large	<118	290

Packaging & Storage	
Packaging	Case 10 boxes of 100
Box Size (mm)	258 x 125 x 82
Carton Dims (mm)	425 x 265 x 273
Shelf Life	Expiration date is 5 years after production date
Storage Instructions	Keep out of direct sunlight, store in a cool & dry place. Keep away from sources of ozone or ignition.
Accreditations	EN ISO 374-1:2016+AI:2018/ Type B; EN ISO 374-5:2016; EN ISO 21420:2020; Compliant with MD regulation (EU) 2017/745 (Class I), EN455-1:2020, EN455-2:2015, EN455-3:2015, EN455-4:2009; PPE Regulation (EU) 2016/425 (CAT III), EN ISO personal protective equipment, as amended to apply in GB.  *Product meets standard EN421 but awaiting official test certification