



Specification

1.32369.0001 Anaerocult® A mini Gas generator system for the incubation of one to four petri dishes in an anaerobic atmosphere

Specification

Growth (Clostridium tetani ATCC 19406)	good to very good
Growth (Clostridium sporogenes ATCC 19404 (WDCM 00008))	good to very good
Growth (Clostridium perfringens ATCC 10543)	good to very good
Growth (Clostridium innocuum ATCC 14501)	good to very good
Growth (Bacteroides fragilis ATCC 25285)	good to very good

Incubation: 24 hours at 35 °C
Reference medium: blood agar

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.